

# Material Safety Data Sheet

PRODUCT CODE NUMBERS:

MDS094186

## SECTION 1

ISSUE DATE: October 21, 2009

IDENTITY: **TOWELETTE, OBSTETRICAL, ALCOHOL FREE**

### MARKETED OR DISTRIBUTED BY:

Medline Industries, Inc.  
One Medline Place  
Mundelein, IL 60060  
**1.800.MEDLINE**

*Emergency Telephone Information:*

**Contact Your Regional Poison Control Center**

## SECTION 2 - HAZARDOUS INGREDIENTS/IDENTITY INFORMATION

Hazardous Components (Specific chemical identity, common name(s))	CAS #	OSHA PEL	ACGIH TLV	Other limits Recommended	% (Optional)
Benzalkonium Chloride		None Established			0.4%

## SECTION 3 - PHYSICAL/CHEMICAL CHARACTERISTICS

Boiling point: 212F	SP Gravity (water=1): 1
Vapor pressure (mm Hg): Not Available	Melting Point: Not Applicable
Vapor density (air=1): Not Available	Evaporation Rate (butyl acetate=1): Not Available
Solubility in water: Soluble	
Appearance and odor: Clear Solution, Citrus Fragrance	

## SECTION 4 - FIRE AND EXPLOSION HAZARD DATA

Flash point (method used): >200F		
Explosive limits: Not Applicable	LEL: Not Applicable	UEL: Not Applicable
Extinguishing media: Water spray, dry chemical, carbon dioxide, alcohol foam		
Special Fire Fighting Procedures: NA		
Unusual Fire and Explosive Hazards: None		

## SECTION 5 - REACTIVITY DATA

Stability: Stable
Conditions to avoid: NA
Incompatibility (materials to avoid): Strong oxidizers, acids and alkalis
Hazardous decomposition or byproducts: Carbon Dioxide, Carbon Monoxide, toxic products of combustion
Hazardous polymerization: Will not occur

# Material Safety Data Sheet

PRODUCT CODE NUMBERS:  
MDS094186

## SECTION 6 - HEALTH HAZARD DATA

Health Hazards (Acute and Chronic): Skin Contact: Can cause irritation of skin Eye Contact: Can cause irritation of the eyes Inhalation: May be irritating to upper respiratory tract.
Medical conditions generally aggravated by exposure: Preexisting skin disorders
<b>Emergency and First Aid procedures:</b> <b>Skin contact:</b> Remove contaminated clothing. Wash affected area with soap and water. Seek medical attention if irritation persists. <b>Eye contact:</b> Flush eyes with copious amounts of water. Seek medical attention if irritation persists. <b>Inhalation:</b> Remove to fresh air. Seek medical attention if difficulty breathing occurs.

## SECTION 7 - SPILL, LEAK, AND WASTE DISPOSAL PROCEDURES

Steps to be taken in case material is released or spilled: Absorb liquid with absorbent material. Wear gloves.
Waste Disposal Method: Follow all local, state and federal regulations.
Precautions to be taken in handling and storing: Store away from incompatible materials.
Other precautions: KEEP OUT OF REACH OF CHILDREN

## SECTION 8 - CONTROL MEASURES

Respiratory Protection (specify type) None normally required
Ventilation: Local exhaust recommended; Mechanical (General) recommended
Protective gloves: Chemical resistant gloves
Eye Protection: Safety glasses/Goggles
Other protective clothing or equipment: Any protective equipment or clothing that will minimize contact with this material.
Work/hygienic practices: Good hygienic practices

## SECTION 9 -- ADDITIONAL INFORMATION

Safety and Hazard Ratings: (Minimal = 0; Slight = 1; Moderate = 2; Severe = 3; Extreme = 4)

Health (Acute-Chronic): 1      Flammability: 0      Reactivity: 0      Environment (Air-Water-Soil): 0

## SECTION: 10 REGULATORY INFORMATION

### USA State Regulations:

## SECTION 11 - ADDITIONAL INFORMATION

*The information provided in this Material Safety Data Sheet has been obtained from sources believed to be reliable. Medline Industries, Inc. provides no warranties, either expressed or implied and assumes no responsibility for the accuracy or completeness of the data contained herein.*

### Document Revision History:

Issued:                      October 21, 2009                      New document  
Revised: